

Optimization and validation of a rapid method for lipid determination in *Nannochloropsis oculata*

Lieselot Balduyck¹, Cedrick Veryser¹, Koenraad Muylaert², Imogen Foubert¹

¹Microbial and Molecular Systems Kulak – Lab Food & Lipids, Kortrijk, Belgium

²Biology Kulak – Lab Aquatic Biology, Kortrijk, Belgium

lieselot.balduyck@kuleuven-kulak.be

INTRODUCTION

Microalgae contain high amounts of (neutral) lipids
= promising source for biodiesel and food applications

Lipid content determination

- Traditionally: gravimetric methods**
 - Advantage: reliable, reproducible
 - Drawbacks: time consuming, labor intensive, high amounts of sample required (not suitable for screening)
 - Alternative: fluorometric method with Nile Red**
 - Advantage: faster
 - Remark: a good validation of this method = correlation between fluorometric and gravimetric method on **independent** samples
→ This has almost not been done!
- Mistakes in other studies:**
- Quantification with external standard curve (e.g. triolein)
Problem: Nile Red interacts differently with “free” lipid globules than with lipids inside a cell wall (De La Hoz Siegler *et al.*, 2012)
 - Different samples taken from the same batch at different times = not independent!

Optimization

- Incubation time** (Figure 2)
Measurement after minimum 90 minutes of incubation at room temperature
- Measuring temperature** (results not shown)
At higher measurement temperatures, fluorescence intensity increases (limit = 45°C), but sample has to adapt to the higher temperature in measuring cell → room temperature has been chosen (25°C)
- NR concentration** (Figure 3)
Fluorescence intensity reaches a maximum at Nile Red concentration of 0,5 mg/mL
- Glycerol addition** (Figure 4)
- Addition of glycerol improves fluorescence intensity
- Maximal fluorescence at 0,25 mL




Figure 2: Fluorescence intensity (excitation wavelength 486 nm, emission wavelength 570 nm) in function of time after addition of Nile Red (at room temperature)




Figure 3: Influence of Nile Red concentration and OD on Nile Red fluorescence intensity (glycerol volume and volume of Nile Red solution remained constant)




Figure 4: Influence of glycerol volume on fluorescence intensity (volume and concentration of NR solution remained constant)

OPTIMIZED PROTOCOL

- 3 mL of algal suspension of OD₇₅₀=0,2
- Spectrofluorometric measurement of blank sample (excitation at 486 nm, emission spectrum from 400 to 800 nm, measurement temperature 25°C)
- Addition of 0,25 mL of glycerol
- Addition of 10 µL of Nile Red solution in acetone (0,5 mg/mL)
- Incubation at room temperature for minimum 90 minutes
- Spectrofluorometric measurement (excitation at 486 nm, emission spectrum from 400 to 800 nm, measurement temperature 25°C)
- Calculation of both relative peak height and relative peak area (corrected for dry weight)

Validation




Figure 5: Correlation curve between total lipid content determined by the optimized Nile Red assay and by gravimetry





Figure 6: Correlation curve between neutral lipid content determined by the optimized Nile Red assay and by gravimetry

- Good correlation between total and neutral lipid content determined by fluorometric and gravimetric methods
- Slightly better correlation with neutral lipid content ($R^2=0,9448$) than with total lipid content ($R^2=0,9144$), but no big difference because total and neutral lipid content are also correlated in *N. oculata* ($R^2=0,9625$)
- Peak height gives the same trends as peak area and comparable correlation coefficients

EXPERIMENTAL SETUP

- Optimization of the method**
 - Incubation time
 - Measuring temperature
 - Nile Red concentration
 - Addition of glycerol: improves penetration of the dye through the cell wall (Doan *et al.*, 2011)
 - Validation: variation of**
 - Nitrogen concentration in medium
 - Harvesting time
- Independent samples



- Total lipid content:** chloroform/methanol extraction according to Ryckebosch *et al.* (2012)
- Neutral lipid content:** Solid Phase Extraction (SPE) into three lipid fractions: neutral lipids, glycolipids and phospholipids

RESULTS AND DISCUSSION

Emission spectra




Figure 1: Emission spectra between wavelengths 400 and 800 nm (excitation wavelength = 486 nm)

Relative peak area = area (*N. oculata* with Nile Red) – area (*N. oculata* blank without Nile Red) – area (Nile Red solution)
Relative peak height = height (*N. oculata* with Nile Red) – height (*N. oculata* blank without Nile Red) – height (Nile Red solution)

CONCLUSION

General conclusions

- Penetration of the dye through the cell wall of *N. oculata* was improved by addition of glycerol in an optimal concentration. This leads to an increased fluorescence intensity.
- A good correlation exists between lipid content determined by the optimized Nile Red assay and by gravimetry ($R^2=0,9448$).
- Gravimetric methods continue to be the most reliable methods for lipid determination, but the Nile Red assay can be used as a screening method (provided that a correlation curve for the used species has been made first!).

Future research

- Optimization and validation for *Isochrysis* is ongoing. The optimized protocol for *N. oculata* can be used, but addition of glycerol has less effect (as *Isochrysis* has no strong cell wall). Validation seems to give another correlation curve.
- Determine if correlation curve is also strain dependent or only species dependent