

## Selective Knockdowns in Maize by Sequence-Specific Protein Aggregation

Camilla Betti, Joost Schymkowitz, Frederic Rousseau,  
and Eugenia Russinova

### Abstract

Protein aggregation is determined by 5–15 amino acids peptides of the target protein sequence, so-called aggregation-prone regions (APRs) that specifically self-associate to form  $\beta$ -structured inclusions. The presence of APRs in a target protein can be predicted by a dedicated algorithm, such as TANGO. Synthetic aggregation-prone proteins are designed by expressing specific APRs fused to a fluorescent carrier for stability and visualization. Previously, the stable expression of these proteins in *Zea mays* (maize) has been demonstrated to induce aggregation of target proteins with specific localization, such as the starch-degrading enzyme  $\alpha$ -glucan water dikinase, giving rise to plants displaying knockdown phenotypes. Here, we describe how to design synthetic aggregation-prone proteins to harness the sequence specificity of APRs to generate aggregation-associated phenotypes in a targeted manner and in different subcellular compartments. This method points toward the application of induced targeted aggregation as a useful tool to knock down protein functions in maize and to generate crops with improved traits.

**Key words** Protein interference, Sequence-specific aggregation, APR, Protein knockdown technology,  $\beta$ -sheet inclusions

---

### 1 Introduction


The global demand for transgenic crops with a high commercial value continuously requires the development of new technologies that are able to circumvent disadvantages or failures of the existing ones. Therefore, a new protein knockdown technology, referred to as protein interference (Pi), initially elaborated in the model species *Arabidopsis thaliana* was further transferred into maize (*Zea mays*) [1]. The Pi technology allows obtaining maize knockdown mutants with a high specificity by inducing a target protein misfolding and aggregation. During the folding process, newly synthesized peptide sequences may take several paths, leading to proper or incorrect protein folding. Misfolded or unfolded proteins can aggregate in insoluble agglomerates [2]. The proposed method is based on the

fact that protein aggregation is mediated by hydrophobic beta-structured interactions in short aggregation-prone segments of polypeptide chains that become exposed upon misfolding, leading to their association into intermolecular aggregates. These short protein sequences are called aggregation-prone regions (APRs) and consist of 5–15 amino acids [3, 4]. APR-mediated protein aggregation is a highly specific self-association process [2, 5].

Importantly, APRs have been shown to be necessary and sufficient for protein aggregation. As most proteins harbor APR segments within their primary structure and because aggregation is sequence specific, it is possible to induce aggregation and subsequently functional depletion of a protein by exposing it to a short target-specific aggregating peptide. The interaction of these APRs depends highly on the sequences and shows a predisposition toward self-association to form intermolecular  $\beta$ -sheets [6]. The presence of APRs within a target peptide sequence can be identified by prediction algorithms, such as TANGO [7].

As most of the *Arabidopsis* and maize proteins harbor aggregation-prone sequence segments within their primary structure, we demonstrated that functional depletion of a protein can be induced by exposing it to a specific APR [1]. Overexpression of different APRs, derived from a single protein or protein family, fused to a fluorescent carrier, results in specific knockdowns similar, albeit often less strong, to previously described genetic mutants. Surprisingly, protein aggregates are not cytotoxic in plants and the synthetic APRs can be targeted also to different cellular compartments or be secreted in the apoplast, allowing highly selective protein knockouts to be obtained. In addition, the intrinsic nature of APRs to form  $\beta$ -sheet structures assures their structural stability when overexpressed in cells. Hence, the APR expression approach can be used as an innovative knockdown method to inactivate proteins by specific in vivo “pull-down” in defined subcellular compartments of maize.

This chapter describes how to select APR peptides for subsequent overexpression in maize to obtain a specific protein knockdown by means of the Pi technology (Fig. 1). An initial scrupulous experimental design is crucial to achieve positive results and several criteria must be followed. The first step is to run the target protein peptide sequence through the TANGO algorithm to predict aggregation-nucleating sequences in target proteins [7]. More than 80% of the *Arabidopsis* proteome contains APRs and, in our experience, similar results can also be expected for the maize proteome [3]. Ideally, the target protein is encoded by a single gene, although also protein families bearing the same APR sequence may be successfully targeted [1]. The APR sequences are selected as the peptide sequences with the highest probability to cross-aggregate by forming intermolecular  $\beta$ -strands. In the TANGO output, this probability is represented as the percentage




**Fig. 1** Protein interference technology working scheme. The target protein primary sequence is processed via TANGO algorithm. Aggregation prone regions (APRs) are selected and engineered in synthetic aggregating blocks (SABs), fused to citrine yellow fluorescent protein (cYFP) carrier. SAB-cYFP stable overexpression in maize leads to the target protein misfolding and functional inactivation

in the  $\beta$ -aggregation propensity calculation. We restricted the in planta analyses to APRs with  $\beta$ -aggregation propensity score above 50%, but lower threshold could also be explored. To ensure the use of APRs that are unique for the target protein, a basic local alignment search tool, BLAST<sup>®</sup>, (<https://blast.ncbi.nlm.nih.gov/Blast.cgi>), is used to compare the APR sequence to all possible fragments of the same size in the maize proteome, provided the word-size is set to be sufficiently small for peptide matching (e.g., [3]).

**Table 1**  
**Nucleotide and peptide sequences for pEN-L1-SAB-L2**

	Nucleotide sequence	Peptide sequence
attL1	CAAATAATGATTTTATTTTGACTGATAGTG ACCTGTTTCGTTGCAACACATTGATGAGCAATGCTTT TTTATAATGCCAACTTTGTACAAAAAAG CAGGCTTA	N/A
attL2	TACCCAGCTTTCTTGACAAAAGTTGGCATT TAAGAAAGCATTGCTTATCAATTTGTTGCAACGAACA GGTCACTATCAGTCAAAATAAAATC ATTATTTG	N/A
Linker	GCCGGCAGCCCGAAGGGCGCCCCGGCCGCCA AGGGCAGCGGCGCC	AGSPKGAPAAKGSGA
N-ter FR	GAGAACGCCGTGGAC	ENAVD
C-ter FR	GGCACCCCGACCGAGGGAGGAG	GTPTREE



**Fig. 2** Schematic representation of constructs expressing different APR variants, indicated as synthetic aggregating blocks (SABs). SABs are fused to cYFP at their C-terminus. *p35S* CaMV 35S promoter, *pEF1α* elongation factor promoter, *N-ter* first 6 AA of the target protein N-terminus, *APR* aggregation prone region, *N-ter FR* APR N-terminal flanking region, *C-ter FR* APR C-terminus flanking region, *cYFP* citrine yellow fluorescent protein, *Kan<sup>R</sup>/Basta<sup>R</sup>* kanamycin/Basta resistance gene, *RB* right border, *LB* left border

Once candidate APRs that are unique for the target protein are selected, their sequence is further optimized by addition of charged flanking regions (FRs) and linkers to obtain synthetic aggregating blocks (SABs) (Table 1). FRs are composed of five to seven charged residues inserted at the N-terminal and C-terminal flanks of the APRs to slow down aggregation and, thus, limit the formation of insoluble aggregates, working as natural “aggregation gatekeepers” [8]. APRs and FRs are cloned in tandem repeats, separated by a linker sequence, to amplify their aggregation potential and are preceded at their 5′ prime by the sequence coding for the first 6 amino acids of the target protein (Fig. 2 and Table 1). Each SAB is also C-terminally attached to a citrine yellow fluorescent protein (cYFP) for

visualization, solubility increase, and biochemical detection. Additionally, it is also possible to express APRs in specific subcellular compartments by fusing a targeting sequence of choice either N-terminally to SABs (for instance, for chloroplast targeting) or C-terminally to cYFP. SABs are finally cloned in Gateway® expression vectors for *Nicotiana benthamiana* and maize transformation and fused to overexpression or tissue-specific promoters [1].

The next step is the evaluation of APR-induced aggregate formation in a transient expression system, i.e., by *Agrobacterium tumefaciens* transformation and infiltration in young *N. benthamiana* leaves. *Agrobacterium*-mediated transient expression is a convenient, fast, and reproducible tool to produce recombinant proteins in approximately 6 days from the bacterial inoculation. This procedure allows confirming, with biochemical assays, the ability of the selected APRs to aggregate. In addition, APR expression, subcellular localization, and tendency to induce aggregates may also be evaluated by confocal laser scanning microscope (CLSM) imaging.

The aggregate formation is biochemically assessed by total protein extraction from *N. benthamiana* under native conditions 3 days after infiltration, followed by Blue Native Polyacrylamide gel electrophoresis (BN-PAGE) and Western blotting with anti-green fluorescent protein (GFP) antibodies. This assay confirms the presence of APR-induced high-molecular weight complexes, indicating that soluble aggregates are formed and will likely target the desired protein of interest.

In BN-PAGE, proteins remain in their native state and, therefore, separate on the basis of their charge-to-mass ratio. Although this technique does not provide direct measurements of the molecular weight, it can be useful to acquire information, such as protein charge or subunit composition, i.e., aggregates formation. Only APRs coded by constructs that induce the formation of complexes with a molecular mass of more than 55–60 kDa (such as the SAB-YFP monomeric form) are selected for subsequent transformation into maize.

The protocol accurately describes the biomolecular procedures required to reach this selection step. Lastly, immature embryos of the maize B104 inbred line [9] are co-cultivated with the *A. tumefaciens* strain EHA101 engineered with the pBb7m34GW expression vectors containing the selected APRs [10]. The derived T0 transgenic maize plants, expressing the APRs at the correct localization, are selected and propagated, ideally until the homozygous state, by routine procedures. The resulting transgenic lines are then analyzed for the appearance of phenotypes compatible with the desired target protein knockdown.

## 2 Materials

**2.1 Plant Material** *Nicotiana benthamiana* plants approximately 4–5 weeks old.

- 2.2 Strains and Competent Cells**
1. *Agrobacterium tumefaciens* C58 strain.
  2. *A. tumefaciens* C58 strain containing the pBin61-P19 vector [11].
  3. *A. tumefaciens* EHA101 strain for maize transformation.
  4. *A. tumefaciens* C58 electro-competent cells.
  5. *Escherichia coli* DH5 $\alpha$  heat shock-competent cells.
  6. EHA101 heat shock-competent cells.

- 2.3 Vectors for Multisite Gateway<sup>®</sup> Cloning**
1. Entry vectors at a concentration of approximately 30–50 ng/ $\mu$ L: pEN-L4-2-R1, containing the cauliflower mosaic virus 35 (CaMV 35S) promoter [12] for overexpression in *N. benthamiana*; pEN-L4-EF1 $\alpha$ -R1 containing the Elongation Factor (EF1 $\alpha$ ) promoter for overexpression in maize (*see Note 1*) [10]; pEN-L1-SAB-L2 vector (custom-made plasmid delivered as a pUC57 vector, *see* Subheading 3.1); pEN-R2-cYFP-L3 to generate translational fusions with cYFP (kind gift of Hilde Nelissen, VIB-Gent University).
  2. Destination vectors at a concentration of approximately 100–150 ng/mL: pK7m34GW for overexpression in *N. benthamiana* and pBb7m34GW for maize transformation [12], conferring in planta resistance to kanamycin and phosphinothricin (BASTA) respectively; these vectors also contain the selectable marker cassette for spectinomycin resistance in bacteria.
  3. Gateway<sup>®</sup> LR Clonase<sup>™</sup> II Plus enzyme mix (Invitrogen).

- 2.4 Media**
1. LB medium: 10 g/L tryptone, 5 g/L yeast extract, 5 g/L NaCl, pH 7.0, autoclave.
  2. Solid LB medium: LB medium supplemented with 15 g/L agar, autoclave.
  3. YEB medium: 5 g/L bacto beef extract, 1 g/L bacto yeast extract, 5 g/L peptone, pH 7.2, autoclave.
  4. Solid YEB medium: YEB medium supplemented with 15 g/L agar, autoclave.

- 2.5 Solutions, Reagents, and Buffers**
1. Antibiotics stock solutions: spectinomycin 100 mg/mL, filter sterilize; gentamicin 20 mg/mL, filter sterilize; rifampicin 50 mg/mL in dimethyl sulfoxide (DMSO); kanamycin 50 mg/mL, filter sterilize; all solutions can be kept at –20 °C in 0.5 mL aliquots (*see Note 2*).

2. TE buffer: 1 mM ethylenediaminetetraacetic acid (EDTA), 10 mM tris(hydroxymethyl)aminomethane (Tris)-HCl, pH 8.0, autoclave.
3. Proteinase K solution.
4. Acetosyringone (3',5'-dimethoxy-4'-hydroxyacetophenone): 0.1 M in DMSO, store at +4 °C.
5. 2-(*N*-morpholino)ethanesulfonic acid (MES) buffer: 0.5 M (pH 5.6), filter sterilize and store at +4 °C (*see* **Note 3**).
6. MgCl<sub>2</sub>: 1 M, autoclave.
7. Infiltration buffer: 10 mM MgCl<sub>2</sub>, 10 mM MES (pH 5.6), 0.1 mM acetosyringone.
8. Protein native extraction buffer (NEB): 20 mM Tris-HCl (pH 7.5), 150 mM NaCl, 18 mM 3-[(3-Cholamidopropyl)dimethylammonio]-1-propanesulfonate hydrate (CHAPS), cComplete™ ULTRA tablets (1 tablet for 10 mL solution).
9. Semi-dry blotting buffer: 12 mM Tris, 96 mM glycine, 20% (v/v) methanol, 0.05% (v/v) sodium dodecyl sulfate (SDS), store at +4 °C.
10. Phosphate-buffered saline (PBS) buffer: 137 mM NaCl, 2.7 mM KCl, 10 mM Na<sub>2</sub>HPO<sub>4</sub>, 1.8 mM KH<sub>2</sub>PO<sub>4</sub>, pH 7.4, autoclave.
11. Washing buffer (PBS-T): 0.1% (v/v) polysorbate20 (Tween20) in PBS.
12. Blocking buffer: PBS-T supplemented with 3% (v/v) skimmed milk.
13. Acetic acid: 8% (v/v).
14. Anti-GFP monoclonal antibody (mouse) horseradish peroxidase (HRP) conjugated (anti-GFP-HRP) (Milenyi Biotec).

## 2.6 Kits, Gels, and Membranes

1. QIAprep Spin Miniprep Kit plasmid purification kit.
2. Quick Start™ Bradford 1× Dye Reagent.
3. NativePAGE™ Novex® 3–12% Bis-Tris Gels 1.0 mm, 10 well (Invitrogen).
4. NativePAGE™ 5% G-250 Sample Additive (Invitrogen).
5. NativePAGE™ running buffer kit (Invitrogen).
6. NativeMark™ unstained protein standard (Invitrogen).
7. Immobilon-P polyvinylidene fluoride (PVDF) membrane: 0.45 µm.
8. Whatman™ 3MM blotting paper.
9. Western Lightning® Plus-Enhanced Chemiluminescence (ECL) substrate.
10. Hyperfilm ECL.
11. 3M™ Micropore™ surgical tape.

## 2.7 Equipment

1. Mortars and pestles.
2. Small (8 cm) and large (13 cm) pots.
3. Vermiculite.
4. Fertilizers: Osmocote<sup>®</sup> and Scotts<sup>®</sup>.
5. Incubators and shakers at 28 °C and 37 °C.
6. Orbital table-top shaker.
7. 1 mL syringes.
8. 1.5 mL microcentrifuge tubes, 50 mL Falcon tubes, 1.5 mL and 2 mL Eppendorf tubes.
9. Protective glasses.
10. Spectrophotometer.
11. Xcell SureLock™ Mini-Cell (Invitrogen) PAGE running system.
12. Semi-dry blotting apparatus.
13. Developing cassettes and film-developing apparatus with developer and fixer.
14. Computer (Windows or Machintosh) system for bioinformatics analyses, databases searches, and sequence analysis.

## 2.8 Databases and Software

1. GDB database (<http://www.maizegdb.org/>).
2. PLAZA2.5 ([http://bioinformatics.psb.ugent.be/plaza/versions/plaza\\_v2\\_5/](http://bioinformatics.psb.ugent.be/plaza/versions/plaza_v2_5/)) [13].
3. TANGO (<http://www.switchlab.org/bioinformatics/tango>) [7].
4. BLASTp<sup>®</sup> (<https://blast.ncbi.nlm.nih.gov/Blast.cgi>).

---

## 3 Methods

### 3.1 APRs In Silico Analysis

1. Select the target protein of interest, preferentially encoded by a single gene (*see Note 4*).
2. Download the peptide sequence of the desired target protein from the maize GDB database or PLAZA2.5 (*see Note 5*).
3. Open TANGO and log in.
4. Select “calculation” in the toolbox, paste the target peptide sequence, and press “Calculate” (*see Note 6*).
5. Keep default parameters and conditions.
6. Select the short APR peptide sequences as the amino acids (AAs) with a cross  $\beta$ -aggregation propensity higher than 50% (*see Note 7*).



7. Run the 5–10 AA APR sequences in BLASTp<sup>®</sup> with the non-redundant (nr) protein sequence database and *Zea mays* (taxid:4577) as organism.
8. Select only APRs that are unique for the target protein to knock down (*see* **Notes 4** and **8**).
9. Design the SABs for LR multisite Gateway<sup>®</sup> cloning in a sequence analysis software (such as Vector NTI<sup>®</sup>, Invitrogen) as represented in Fig. 2. Insert the following synthetic sequences in this order: six AAs of the N-terminus (N-ter) of the target protein; five AAs of N-ter flanking region (N-ter FR); APR; seven AAs of the C-terminus (C-ter) FR; Linker; repeat the sequence, N-ter FR-APR-C-ter FR-Linker (*see* Fig. 2 and Table 1 for SAB schematic representation and sequences) (*see* **Note 9**).
10. Place an order for the SAB sequences to a gene synthesis company, such as Genescript<sup>®</sup>, after flanking them with the attL1 and attL2 cloning sites, so they are ready-to-use for three fragment recombination multisite Gateway<sup>®</sup> cloning (*see* **Note 10**).
11. Choose the SABs to be delivered in a cloning vector, such as a pUC57 plasmid, for direct use as entry vectors for LR recombination (pEN-L1-SAB-L2).

### **3.2 APRs Cloning in pK7m34GW and pBb7m34GW Expression Vectors**

1. For the MultiSite Gateway<sup>®</sup> LR recombination reaction between multiple entry clones, in a 1.5 mL microcentrifuge tube, mix 6  $\mu$ L TE, 1  $\mu$ L pK7m34GW or pBb7m34GW (approximately 100–150 ng/ $\mu$ L), 1  $\mu$ L pEN-L4-2-R1 or 1  $\mu$ L pEN-L4-EF1 $\alpha$ -R1, 1  $\mu$ L pEN-L1-SAB-L2 and 1  $\mu$ L pEN-R2-cYFP-L3 (*see* **Note 11**) and keep at room temperature.
2. Briefly vortex the LR Clonase<sup>™</sup> II Plus enzyme mix and keep it on ice.
3. Add 2  $\mu$ L of LR Clonase<sup>™</sup> II Plus enzyme mix to the microcentrifuge tube (**step 1**) and gently mix.
4. Spin down and incubate the LR reaction at 25 °C for 16 h (or overnight).
5. Add 1  $\mu$ L of 2  $\mu$ g/ $\mu$ L Proteinase K solution to neutralize the clonase enzyme.
6. Gently mix and spin down briefly.
7. Incubate the reaction at 37 °C for 10 min, spin down, and place on ice.
8. Transform 2  $\mu$ L of the LR reaction mix into 50  $\mu$ L *E. coli* DH5 $\alpha$  heat shock-competent cells.

9. Spread 20  $\mu$ L and 100  $\mu$ L aliquots onto plates of solid LB medium containing 100  $\mu$ g/mL spectinomycin.
10. Incubate the plates upside down overnight at 37 °C.
11. Select six single *E. coli* colonies and inoculate them in 10 mL of LB medium containing 100  $\mu$ g/mL spectinomycin.
12. Incubate overnight at 37 °C in a shaker at 220 rpm.
13. Isolate the expression vectors with the QIA miniprep kit, including all clean-up steps.
14. Elute the DNA with 50  $\mu$ L of sterile bi-distilled water.
15. Determine the DNA concentration with the spectrophotometer.
16. Sequence six isolated plasmids per construct.
17. After the sequence analysis confirmation, make glycerol stocks and transform one plasmid per construct into *Agrobacterium* C58 electrocompetent cells or EHA101 heat shock-competent cells.
18. Plate the transformed bacteria onto YEB agar plates containing 20  $\mu$ g/mL gentamycin, 100  $\mu$ g/mL spectinomycin, and 25  $\mu$ g/mL rifampicin for C58 *Agrobacterium* strain, or onto YEB agar plates containing 100  $\mu$ g/mL spectinomycin and 25  $\mu$ g/mL kanamycin for the EHA101 *Agrobacterium* strain.
19. Incubate upside down at 28 °C for 2–3 days (in the dark for plates containing rifampicin).
20. Pick single colonies, streak them on selective YEB agar, and make glycerol stocks (*see* **Notes 12** and **13**).

### 3.3 APRs Biochemical Assessment by BN-PAGE and Western Blotting

#### 3.3.1 *Agrobacterium*- mediated Infiltration of *p35S-SAB-cYFP Constructs*

1. Sow approximately 20–30 *N. benthamiana* seeds in a small pot containing wet soil and transfer for germination in a greenhouse under controlled growth conditions at 21 °C and 65% humidity.
2. After 2 weeks, transfer seedlings in single, large pots containing soil and vermiculite in a 3:1 ratio mixed with a teaspoon full of Osmocote® fertilizer.
3. Water and keep at 21 °C and 65% humidity under controlled growth conditions.
4. Fertilize after 2 weeks with 1 g/L Scotts fertilizer 20:20:20.
5. Inoculate *Agrobacterium* strains transformed with the p35S-SAB-cYFP constructs, in 10 mL YEB supplemented with 20  $\mu$ g/mL gentamycin, 100  $\mu$ g/mL spectinomycin, 10 mM MES (pH 5.6), and 20  $\mu$ M acetosyringone, in 50-mL Falcon® tubes leaving the cap half open and sealed with 3M™ Micro-pore™ surgical tape (*see* **Note 2**).

6. In parallel, inoculate *Agrobacterium* containing the pBin61-P19 vector in 10 mL YEB added with 20 µg/mL gentamycin, 50 µg/mL kanamycin, 10 mM MES (pH 5.6), and 20 µM acetosyringone (*see Note 2*). Coexpression of p19 protein from tomato bushy stunt virus is used for suppression of transgene silencing resulting in increased expression of the desired proteins [11].
7. Incubate the cultures at 28 °C for 2 days.
8. To measure and record the optical density at 600 nm (OD<sub>600</sub>) of the cultures, make a 1:5 dilution of each culture and place 1 mL of the diluted cultures into a plastic cuvette. Read the OD with the spectrophotometer. Make blank controls with YEB medium and appropriate antibiotics.
9. Calculate the amount of each culture needed to have a final OD<sub>600</sub> of 1.5 in a final volume of 2 mL (*see Note 14*).
10. Transfer the bacterial cultures in 2 mL Eppendorf tubes and spin down the appropriate volumes by centrifugation at 4000 × *g* for 20 min.
11. Discard the supernatant.
12. Resuspend the pellets obtained in **step 10** in 2 mL infiltration buffer by vigorous flicking (*see Note 3*).
13. Incubate at room temperature for 2–3 h to let the acetosyringone activate the *Agrobacterium* virulence genes.
14. After incubation, mix the *Agrobacterium* cultures 1:1 and make a 1 mL infiltration mix (*see Note 15*).
15. Prepare to infiltrate the bacterial mix into the abaxial side of approximately 4-week-old *N. benthamiana* leaves with a 1 mL syringe without needle. Make a small incision in the epidermis on the leaf abaxial side by a soft touch with a hypodermic needle.
16. To infiltrate the *Agrobacterium* cell mixture, simply press the syringe on the leaf at the incision level and exert a counter-pressure with your finger tip on the other side (*see Note 16*).
17. Highlight the perimeter of the injected area that will turn dark green with a permanent marker on the upper side of the leaf.
18. Continue to inject until the infiltration mixture will stop to expand.
19. Repeat the procedure until the entire leaf is infiltrated.
20. Transfer the plants under controlled growth conditions at 21 °C and 65% humidity and water them.
21. Three days after infiltration, collect the injected leaf material, flash freeze in liquid nitrogen, and store at –70 °C for subsequent biochemical analyses (*see Note 17*).

## 3.3.2 BN-PAGE

1. Grind the leaf material obtained in Subheading 3.3.1 in ice-cold mortars to a fine powder.
2. Fill an ice-cold 2 mL Eppendorf tube containing 500  $\mu$ L NEB with the ground material.
3. Vortex and place on ice.
4. Incubate approximately 15–20 min on ice and after 10 min mix thoroughly by inverting the tubes.
5. Centrifuge at  $14,000 \times g$  for 20 min at 4 °C to pellet the cell debris.
6. Collect the supernatant in a new ice-cold 1.5 mL Eppendorf tube and keep on ice. If leaf particulate is still present in the supernatant, repeat centrifugation in **step 5** for 10 min (*see Note 18*).
7. To measure protein concentration, perform a Bradford assay [14] by preparing an appropriate number of cuvettes containing 990  $\mu$ L of Quick Start™ Bradford 1 $\times$  Dye Reagent (*see Note 19*).
8. Add to the cuvette 10  $\mu$ L of protein extracts obtained in **step 6**, diluted 1:10 in NEB and add 10  $\mu$ L of NEB in the blank control.
9. Mix thoroughly and incubate 15 min on the work-bench.
10. After incubation, record the OD at 595 nm (OD<sub>595</sub>), after putting the blank control at a null value.
11. Calculate the protein concentrations by using the calibration curve extrapolated from bovine serum albumin (BSA) OD<sub>595</sub> reads.
12. For the BN-PAGE electrophoresis, prepare 1.5 mL Eppendorf tubes containing the samples obtained in **step 6** diluted in NEB to have 60  $\mu$ g of total proteins in a 25  $\mu$ L volume.
13. Add 10% (v/v) glycerol and 2  $\mu$ L NativePAGE™ Coomassie G-250 Sample Additive to each sample (*see Note 20*).
14. Wash the wells of a 3–12% Novex® Bis-Tris gradient gel with NativePAGE™ cathode buffer twice prior to assembling into the BN-PAGE X-cell Sure-lock™ running apparatus (*see Note 20*).
15. Prepare the Xcell SureLock™ Mini-Cell running apparatus: mount the gels, fill the inner chamber with 200 mL NativePAGE™ cathode buffer (10 mL NativePAGE™ 20 $\times$  running buffer, 1 mL NativePAGE™ 20 $\times$  cathode buffer additive in sterile water) and the outer chamber with 400 mL NativePAGE™ anode buffer obtained by diluting 20 $\times$  NativePAGE™ running buffer.

16. Load the samples prepared in **step 12** into a gradient Native-PAGE™ Novex® 3–12% Bis-Tris gel.
17. Use 10 µL of NativeMark™ unstained protein standard as a marker, added with 10% (v/v) glycerol and 2 µL of Native-PAGE 5% G-250 sample additive.
18. Perform the electrophoresis at room temperature for 2 h at 120V (*see Note 21*).

### 3.3.3 Western Blotting

1. Activate a PVDF membrane, previously cut to the same size of the gel, by incubation in 100% (v/v) methanol for 15 s.
2. Equilibrate the membrane in semi-dry blotting buffer for 5 min.
3. Soak five Whatman™ filter papers, previously cut to the same size of the gel, and the PVDF membrane in semi-dry blotting buffer.
4. Transfer the proteins in the gel from Subheading 3.3.2 onto a PVDF membrane by semi-dry blotting.
5. After the run, carefully dismantle the gel from the plastic cast and place it in a semi-dry blotting apparatus with material in **step 1** in this order: from anode to cathode, two Whatman™ filter papers, PVDF membrane, native gel, and three Whatman™ filter papers.
6. Perform the blot at a fixed current of 45 mA for 1.15 h (*see Note 22*).
7. Fix proteins on the PVDF membrane in 8% (v/v) acetic acid for 20 min on a table-top shaker.
8. Air-dry the membrane.
9. Destain the Coomassie blue from the PVDF membrane with 100% (v/v) methanol and wash twice in PBS buffer to remove the methanol residues.
10. Block the proteins on the membrane by incubation in blocking buffer overnight at 4 °C or for 1 h at room temperature (*see Note 23*).
11. Discard the blocking solution and add 10 mL of blocking solution containing 1:5000 anti-GFP-HRP monoclonal antibody for immunodetection of the SAB-cYFP aggregates (*see Note 16*).
12. Incubate the membrane for 1 h at room temperature.
13. Discard the antibody solution (*see Note 24*) and wash the membrane with PBS-T buffer for 15 min.
14. Discard the solution and wash two times with PBS-T for 5 min.
15. Gently dry the membrane for 2 s with adsorbent paper to remove excess of wash solution.

16. To develop the membrane with an ECL kit, mix equal amounts of enhanced luminol reagent and oxidizing reagent (usually 1 mL of each substrate is sufficient) and incubate for 1 min at room temperature.
17. Gently dry the membrane with adsorbent paper for 2 s and place it in a developing cassette between transparent plastic foils.
18. In a dark room, place a Hyperfilm ECL on the top of the membrane with firm hands.
19. Detect the chemiluminescence emitted from the filter after 5 min by using a film developer apparatus.
20. Repeat the procedure at different times, usually between 5 and 60 min depending on the cYFP signal intensity.
21. Analyze the results and select the APRs that induce high molecular-weight aggregates, when expressed in the 35S-SAB-YFP constructs, for subsequent maize transformation (*see* **Note 25**).

---

## 4 Notes

1. To drive transgene overexpression in maize, better results are obtained when the elongation factor promoter (pEF1 $\alpha$ ) from *Brachypodium distachyon* is used instead of the p35S from CaMV [10]. When a mesophyll-specific expression is desired, a phosphoenol pyruvate carboxylase (PEPC) promoter can be successfully used [1].
2. Rifampicin is a light-sensitive antibiotic that degrades fairly quickly in the presence of light and, in solution, it loses its activity. Therefore, because shakers in most laboratories are kept under light conditions, we recommend avoiding addition of this antibiotic to the YEB medium. The *Agrobacterium* strain transformed with the pBIN61-P19 vector is resistant to rifampicin, gentamycin, and kanamycin, but, for the above-mentioned reasons, only the last two are used for inoculation in liquid YEB medium. Moreover, when preparing the medium for P19 cultures, consider an amount that will be sufficient to be injected together with each p35S-SAB-cYFP *Agrobacterium*-expressing construct. Finally, in the experiment, always include also a negative culture control (medium with antibiotics) to make sure to avoid any medium contamination.
3. The pH of the MES solution can be unstable. If the solution is not freshly prepared, always control that the pH is in the 5.6–5.7 range before use. Filter sterilization is advised, because after autoclaving the solution can turn yellowish and, although

the pH does not change measurably, the nature of the by-products is unknown.

4. This approach could be applied also to knock down protein families, when predicted aggregation-prone peptides would overlap across all family members [1]. To be able to evaluate the results in the optimal way possible, it would be ideal to knock down a protein family of which the knockdown would lead to known, or easy-to score, phenotypes.
5. The PLAZA2.5 bioinformatics platform is ideal when searching for maize protein sequences orthologous to others from different crop or model plant species [13].
6. In the calculation output, it is represented in the “ $\beta$  aggregation propensity” column.
7. The TANGO algorithm available online is suitable only to process proteins with a sequence of up to 500 amino acids. When longer sequences need to be analyzed, please contact the SWITCH laboratory (VIB-KU Leuven, Belgium).
8. Our results show that APRs with similar TANGO scores ( $\geq 50\%$ ) and one mismatch, i.e., a sequence identity around 90%, can bind the target proteins, whereas APRs with a TANGO score  $\leq 50\%$ , but bearing two mismatches, i.e., less than 80% sequence identity, do not interact. Therefore, it is advised to select APRs with a TANGO score  $> 50\%$  and unique for the target protein, unless a downregulation of an entire protein family sharing a similar APR sequence is desired [1].
9. Nucleotide sequences must be optimized for expression in maize.
10. Before placing the order for the SAB sequences, make sure that the attL1 e attL2 recombination sequences for the three fragment multisite Gateway® cloning are correct. In fact, we noticed that when the Gateway® entry clones are recombined in silico, the attL1 and attL2 nucleotide sequences contain one or two nucleotides that are not considered recombination sequences by the software, and/or are not listed in the manual. When these nucleotides are missed in the synthetic sequences, it will obviously lead to incorrect translation of the APRs due to a frameshift. Therefore, we advise to include as attL1 and attL2 sequences all the nucleotides inserted after the in silico recombination between the 3' end of the promoter sequence and the 5' of the SAB. As the SAB sequences are highly hydrophobic and, hence, problematic to synthesize, it is recommended to use a well-known and skilled gene synthesis company, such as Genscript®. After the delivery of the order, follow the instructions provided with the supplied plasmids to dissolve the DNA at a working concentration of usually 200 ng/ $\mu$ L.

11. The concentrations of all entry vectors used are at approximately 30–50 ng/ $\mu$ L.
12. It is possible to further confirm the *Agrobacterium* clones by colony PCR or miniprep followed by sequencing, although, in our experience, these steps can be quite troublesome to carry out on *Agrobacterium* due to the thick cell walls. For this reason, the initial steps of the cell lysis must be extended and improved in both the cases.
13. Recombinant C58 *Agrobacterium* containing the p35S-SAB-cYFP in pK7m34GW expression vectors will be used for *Agrobacterium*-mediated infiltration experiments. Recombinant EHA101 *Agrobacterium* expressing pBb7m34GW vectors containing APRs selected at the end of Subheading 3.3.3, fused to cYFP and under the pEF1 $\alpha$  promoter (or tissue-specific promoters), will be further inserted into immature B104 maize embryos.
14. For injecting two strains, one containing p35S-SAB-cYFP and one the pBin61-P19 vectors, each strain should be at OD<sub>600</sub> of 0.75. It is important that all constructs have the same OD<sub>600</sub>  $\pm$  0.05.
15. Prepare at least 6 mL of infiltration mix when injecting leaves for subsequent biochemical assays. To have enough leaf material for protein extraction, at least one entire leaf per construct combination is necessary; 3 mL infiltration mix should suffice for one leaf injection. We recommend injecting two leaves from different plants with each construct combination to avoid possible differences in expression levels after in planta transformation.
16. When punching the leaves, caution should be taken not to perforate the leaf through both sides; otherwise, the *Agrobacterium* mixture will not infiltrate into the intercellular space, but would just pass through the wound to the other leaf side. When multiple constructs are injected at the same time, be extra careful to not cross-contaminate the samples. Gloves should be changed after each construct combination injection and protective glasses must be worn during all the steps.
17. Although it is not strictly required to carry out the Pi technology main procedures, at this stage, it is also possible to check the APR expression by cYFP CLSM imaging. Although YFP is clearly visible by imaging with GFP settings, for visualization under the best conditions, it is advised to adjust the confocal microscope settings as much as, and whenever, possible for the YFP excitation (peak at 514 nm) and emission spectra (peak at 527 nm). To this aim, the best results, especially when imaging maize leaves, can be obtained by using a confocal microscope with adjustable BP and LP filters (such as the Olympus<sup>®</sup>



FluoView™ FV1000 inverted confocal microscope). Importantly, cYFP is also the best fluorescent probe for maize CLSM imaging in our hands. As YFP protein only carries three amino acid differences with GFP, its use does not affect the extraordinary number of possibilities to perform biochemical analyses by means of commercial anti-GFP antibodies, beads for immunoprecipitation, etc.

18. Flash-frozen tissues can be stored for a long time at  $-70^{\circ}\text{C}$  before protein extraction. However, once total proteins are isolated under native conditions, they need to be processed immediately by BN-PAGE to avoid unwanted denaturation and/or degradation of the protein aggregates contained in the extracts. To this end, samples must be kept also on ice during the whole procedure. The BN-PAGE technique may be also applied to maize leaf material by following the same procedure as described here.
19. Samples must be read in triplicates and a calibration curve with 0–1000  $\mu\text{g}/\text{mL}$  BSA must be done in parallel.
20. In the procedure described here, a NativePAGE™ NOVEX Bis-Tris Gel System (Invitrogen) is used. The NativePAGE™ procedure is based on the use of Coomassie Brilliant Blue G-250 in the sample buffer (sample additive), because it binds to proteins, conferring them a negative charge without denaturation. Moreover, Tricine (negatively charged) and Bis Tris (positively charged) salts are used as electrophoresis ions, allowing the electrophoresis procedure to proceed at nearly neutral pH (7.5–7.7). As the cathode buffer contains Coomassie blue, it has a deep blue color making gel loading very tricky. Therefore, to get a visual help during loading, it is advised to mark the wells in the precast gels with a black marker before starting the Xcell SureLock™ running cassette assembly. Due to the delicate nature of the 3–12% gradient native gels, it is also recommended to use always commercial precast gels together with reagents from the same company to avoid unwanted running problems due to unequal ions distribution, which are very common when gradient gels are cast in the laboratory.
21. Occasionally, the run may stop due to low current. In that case, set the voltage to 180 V for 40 min.
22. Always perform the semi-dry blotting at a low current; otherwise, the native gel could melt.
23. This step, as well as **steps 5–7**, is performed on an orbital table-top shaker.
24. The YFP fluorescent protein is recognized by a great variety of commercial polyclonal and monoclonal anti-GFP antibodies, because it has a primary, secondary, and tertiary structure

extremely similar to that of GFP. We advise to use the anti-GFP HRP-conjugated monoclonal antibody, because it is very specific and also allows the secondary antibody incubation step to be skipped during the immunoblotting experimental procedures. After use, the antibody solution can be kept at  $-20^{\circ}\text{C}$  and reused several times (tested up to four times).

25. After blotting, YFP-labeled aggregates may appear either as definite bands or as smears, at least twice the size of the APR-YFP monomeric form (approximately 55–60 kDa). The smears indicate that differently sized aggregates are formed. In our experience, APRs inducing these kinds of aggregates are as efficient as those that induce only aggregates of one kind in achieving the target protein knockdown. Lastly, when an antibody against the target protein is available, we propose to perform a co-immunoprecipitation procedure, after co-expression of the target protein and the selected APRs in *N. benthamiana*, to confirm that the APRs of choice interact with the target protein of interest. This will strongly support a positive outcome of the Pi knockdown approach.

---

## Acknowledgments

This work was supported by grants from the Agency for Innovation by Science and Technology (“Strategisch Basisonderzoek” project no. 60839), Ghent University (“Industrieel Onderzoeksfonds” F2011/IOF-Advanced121 and F2014/IOF-StarTT261 and Multidisciplinary Research Partnership “Biotechnology for a Sustainable Economy” no. 01MRB510W), the Interuniversity Attraction Poles Program (IUAP VII/29), initiated by the Belgian State, Science Policy Office, University of Leuven, and the European Research Council under the European Union’s Horizon 2020 Framework Programme (ERC Grant agreement 647458). We thank Martine De Cock for proofreading the chapter.

## References

1. Betti C, Vanhoutte I, Coutuer S, De Rycke R, Mishev K, Vuylsteke M, Aesaert S, Rombaut D, Gallardo R, De Smet F, Xu J, Van Lijsebettens M, Van Breusegem F, Inzé D, Rousseau F, Schymkowitz J, Russinova E (2016) Sequence-Specific protein aggregation generates defined protein knockdowns in plants. *Plant Physiol* 171(2):773–787
2. Tyedmers J, Mogk A, Bukau B (2010) Cellular strategies for controlling protein aggregation. *Nat Rev Mol Cell Biol* 11:777–788
3. Rousseau F, Serrano L, Schymkowitz JWH (2006) How evolutionary pressure against protein aggregation shaped chaperone specificity. *J Mol Biol* 355:1037–1047
4. Goldschmidt L, Teng PK, Riek R, Eisenberg D (2010) Identifying the amyloids, proteins capable of forming amyloid like fibrils. *Proc Natl Acad Sci U S A* 107:3487–3492
5. Mitraki A (2010) Protein aggregation: from inclusion bodies to amyloid and biomaterials. *Adv Protein Chem Struct Biol* 79:89–125

6. Ganesan A, Debulpaep M, Wilkinson H, Van Durme J, De Baets G, Jonckheere W, Ramakers M, Ivarsson Y, Zimmermann P, Van Eldere J, Schymkowitz J, Rousseau F (2014) Selectivity of aggregation-determining interactions. *J Mol Biol* 427:236–247
7. Fernandez Escamilla A-M, Rousseau F, Schymkowitz J, Serrano L (2004) Prediction of sequence dependent and mutational effects on the aggregation of peptides and proteins. *Nat Biotechnol* 22:1302–1306
8. De Baets G, Van Durme J, Rousseau F, Schymkowitz J (2014) A genome-wide sequence-structure analysis suggests aggregation gatekeepers constitute an evolutionary constrained functional class. *J Mol Biol* 426:2405–2412
9. Hallauer AR, Lamkey KR, White PR (1997) Registration of five inbred lines of maize: B102, B103, B104, B105, and B106. *Crop Sci* 37:1405–1406
10. Coussens G, Aesaert S, Verelst W, Demeulenaere M, De Buck S, Njuguna E, Inzé D, Van Lijsebettens M (2012) *Brachypodium distachyon* promoters as efficient building blocks for transgenic research in maize. *J Exp Bot* 63:4263–4273
11. Voinnet O, Lederer C, Baulcombe DC (2000) A viral movement protein prevents spread of the gene silencing signal in *Nicotiana benthamiana*. *Cell* 103:157–167
12. Karimi M, Bleys A, Vanderhaeghen R, Hilson P (2007) Building blocks for plant gene assembly. *Plant Physiol* 145:1183–1191
13. Proost S, Van Bel M, Sterck L, Billiau K, Van Parys T, Van de Peer Y, Vandepoele K (2009) PLAZA: a comparative genomics resource to study gene and genome evolution in plants. *Plant Cell* 21:3718–3731
14. Bradford MM (1976) A rapid and sensitive method for the quantitation of microgram quantities of protein utilizing the principle of protein-dye binding. *Anal Biochem* 72:248–254